


Space Saver Header Bars

Space saver header bars are ideal for manifold applications where space is limited. The modular header bars feature dual ports, allowing for optimum cylinder placement. The complete kits include the following components for each side of the manifold: Main shutoff valve; header bar segments with check valves at each connection port; 36" flexible pigtails, except for high pressure oxygen, which are 24" copper; mounting brackets.

The kits can be used for low- or high-pressure gases, except for low pressure oxygen. The header bars can be expanded in the future if additional capacity is required.



Ordering Information

Kit	Oxygen	Nitrous Oxide	Nitrogen	Air	Carbon Dioxide	Length*
1 X 1	261694-01	261695-01	261696-01	261697-01	261698-01	43.3"
2 X 2	261694-02	261695-02	261696-02	261697-02	261698-02	52.1"
3 X 3	261694-03	261695-03	261696-03	261697-03	261698-03	63.1"
4 X 4	261694-04	261695-04	261696-04	261697-04	261698-04	51.7"
5 X 5	261694-05	261695-05	261696-05	261697-05	261698-05	63.1"
6 X 6	261694-06	261695-06	261696-06	261697-06	261698-06	71.9"
7 X 7	261694-07	261695-07	261696-07	261697-07	261698-07	73.1"
8 X 8	261694-08	261695-08	261696-08	261697-08	261698-08	71.4"
9 X 9	261694-09	261695-09	261696-09	261697-09	261698-09	92.9"
10 X 10	261694-10	261695-10	261696-10	261697-10	261698-10	101.7"
11 X 11	261694-11	261695-11	261696-11	261697-11	261698-11	112.6"
12 X 12	261694-12	261695-12	261696-12	261697-12	261698-12	111.4"
13 X 13	261694-13	261695-13	261696-13	261697-13	261698-13	122.4"
14 X 14	261694-14	261695-14	261696-14	261697-14	261698-14	131.2"

* Length is for both sides of header bar connected to an 18" wide manifold cabinet. Total horizontal space will be dependant upon cylinder placement.

Header Bar Segments with Dual Connections

Description	Oxygen	Nitrous Oxide	Nitrogen	Air	Carbon Dioxide	Length
Shutoff valve and 1 outlet	263126-1	263129-1	263128-1	263127-1	263130-1	13.5"
Shutoff valve and 2 outlets	263126-2	263129-2	263128-2	263127-2	263130-2	17.9"
Shutoff valve and 4 outlets	263126-4	263129-4	263128-4	263127-4	263130-4	17.7"
Extension with 2 outlets	263131-2	263134-2	263133-2	263132-2	263145-2	11.5"
Extension with 4 outlets	263131-4	263134-4	263133-4	263132-4	263145-4	11.5"
Extension with 6 outlets	263131-6	263134-6	263133-6	263132-6	263145-6	16.5"
90° Elbow for header bar	261677					

Header Bar Segments for Low Pressure Oxygen

Description	Oxygen	Length
Shutoff valve and 1 outlet	263962	13.5"
Shutoff valve and 2 outlets	263963	17.9"
Shutoff valve and 3 outlets	263964	23"
Extension with 2 outlets	263965	11.5"
Extension with 3 outlets	263966	16.5"

© 2023 Ohio Medical LLC. This document contains information that is proprietary and confidential to Ohio Medical. Use of this information is under license from Ohio Medical. Any use other than that authorized by Ohio Medical is prohibited. Ohio Medical and the Ohio Medical logo are registered trademarks of Ohio Medical. Amvex and the Amvex logo and Aeros are registered trademarks and Squire-Cogswell is a trademark of Ohio Medical. NFPA is a registered trademark of National Fire Protection Association Inc. ISO is a registered trademark of International Organization for Standardization.